Creatine Kinase-MB Test Kit

(Quantitative Fluorescent Immuno-chromatography)

【Product Name】

Generic Name: Creatine kinase-MB Test Kit (Quantitative Fluorescent Immuno-chromatography)

【Package Spec】

25 PCS/box

【Intended Use】

The creatine kinase-MB test kit (Quantitative Fluorescent Immuno-chromatography) is used to quantitatively detect the concentrations of Creatine kinase-MB (CK-MB) in human serum, plasma and whole blood samples in vitro by using fluorescent probe labeled immuno-chromatography technology combined with Immunity analyzer.

Creatine kinase catalyzes the transfer of high-energy phosphate bonds from ATP to creatine to produce creatine phosphate. Creatine kinase is a dimeric enzyme which including four isoenzyme: skeletal muscle CK-MM, brain tissue CK-BB, cardiac CK-MB, and mitochondrial creatine kinase. CK-MB has a molecular weight of 80,000 and consists of two subunits (MW = 40000 per subunit): the subunit M represents muscle type and the subunit B represents brain type. CK-MB is mainly found in myocardial tissue, accounting for 20% of total CK activity, and also has a small amount of distribution in skeletal muscle. Serum CK-MB content is important for the diagnosis of myocardial infarction such as acute myocardial infarction and myocarditis. CK-MB can be detected in circulating blood in the early stage of myocardial injury, and rises 2-6 hours after myocardial injury, peaks within 12-24 hours, then falls to normal level (36-72 hours). The elevation of CK-MB concentration is accompanied by ECG evolution and a history of chest pain. The change in CK-MB content and time is important for assessing the onset time of myocardial injury, infarct volume, and post-reperfusion injury. In addition, skeletal muscle damage, rhabdomyolysis, and stroke could lead to an increase in CK-MB as well.

Clinically, CK-MB is often detected in combination with myoglobin and troponin I. In the diagnosis of myocardial injury-related diseases, troponin I is highly specific, and myoglobin and CK-MB have high early sensitivity. The combined detection of the three can improve the specificity and sensitivity of diagnosis.

【Principle of Testing】

The Test Kit adopts quantum dot fluorescence labeling technology. Mouse anti-human CK-MB monoclonal antibody I was coated on nitrocellulose membrane. Fluorescence labeled mouse anti-human CK-MB monoclonal antibody II was distributed in glass cellulose membrane. If the blood samples contain CK-MB, the analyte can be captured by quantum dot labeled mouse anti-human CK-MB monoclonal antibody II. The conjugate can move to the test line area on the nitrocellulose membrane by capillary action and be captured by the corresponding antibody I contained therein to form a quantum dot fluorescent test line. Within a certain range, the concentration of CK-MB in the sample has a positive correlation with the fluorescence intensity of the test line. The concentrations of CK-MB in the sample can be calculated and determined by the immunofluorescence analyzer.

【Main components】

The main components of the kit include Test Card, Calibration QR Code (on the box), Diluent and Instruction Manual.

The main components of the kit include test card, QR Code Calibration Card, Diluent and Instruction Manual.

Details are as follows:

1. The test card is composed of a test strip and a plastic clip. The main components of the strip are:

- 1) Mouse anti-human Creatine Kinase-MB antibody (immobilized on the "Test Line" area of nitrocellulose membrane);
- 2) Fluorescent labeled mouse anti-human Creatine Kinase-MB monoclonal antibody (immobilized on the glass cellulose membrane);
- 3) Test strip back supporting material (Polyvinyl Chloride, PVC).

2. Calibration QR code: QR Code located on the box is used to carry the product parameters. Please scan the QR code on the box with Immunofluorescence Analyzer for the input of test information.

3. Diluent

Sample diluent solution : contains phosphate buffer and Proclin-300.

Note: Please do not exchange the components in different batches of kits.

【Storage conditions and Period of Validity】

The period of validity of the test kit is 18 months while the package is sealed and stored at 2-30 °C. Once the package is opened, please place the test kit at 20-30 °C. The opened test kit should be used within 30min if the humidity is no less than 20% and below 80%RH.

【Applicable instruments】

The Immunofluorescence Analyzer manufactured by Deniz Biyokimya Ürünleri A.S. $R_{2D_{2}}=V_{2}$

【The requirements of sample】

This kit is suitable for serum/plasma/whole blood samples.

Serum/plasma collection: the serum or plasma should be separated for analysis immediately after blood collection to avoid hemolysis. Use within 24 hours if the sample is stored at room temperature; use within 72 hours if the sample is stored at 2-8 °C; use within 3 months if the sample is stored at -20 °C with one cycle of freeze and thaw.

The collection of whole blood: collect the blood using commercially available blood collection tubes containing EDTA, sodium heparin or sodium citrate anticoagulants. Once the blood is collected, gently invert tube several times to avoid any hemolysis. If the sample cannot be tested within 24 hours, it shall be stored at 2-8°C. Samples being stored more than 3 days should not be used for examination.

Whole blood/plasma is recommended to be collected using the blood collection tubes containing anticoagulant of EDTA (concentration 3.7-5.4 μmol/L), sodium heparin (concentration of 12 IU-30 IU per milliliter of blood) or sodium citrate (concentration 105-109 mmol/L). To separate plasma, it is recommended to centrifuge the blood sample with 500g (3000 rpm) for 5 minutes.

Fresh samples are recommended. Blood samples of hyperlipidemia, jaundice and hemolysis are not recommended.

The sample must reach the equilibrium status at room temperature (20 ~ 30 °C) before use.

【Test Method】

Please read the user manual for test kit and Immunofluorescence Analyzer carefully before use, following the requirements and steps strictly. All the test shall be operated under room temperature (20 ~ 30 °C) and the relative humidity should be no less than 20%RH and below 80%RH.

The signal ratio (T/C) of Test Line (T) and Quality Control Line (C) reading from the immunofluorescent analyzer is converted to analyte concentration using the standard curve information carried by QR code on the box.

The Quality Control Line is the quality control standard of the test kit. If the Quality Control Line has no signal or abnormal signal (exceed the quality control scope of the analyzer), the test result will be invalid, and it is recommended to repeat the test.

Immunofluorescence Analyzer has two detection modes: "quick test" and "standard test". Please select the right detection mode according to user manual of Immunofluorescence Analyzer before detection.

The detail operation step is as belows:

- 1) Allow both the test kits and samples balance to room temperature (20 ~ 30 °C) before performing the test. Confirm that the Immunofluorescence Analyzer is connected to powered supply and turned on before use.
- 2) Before using each lot of the test kit, please select "quick test" mode of the Immunofluorescence Analyzer and QR code on the box to record the lot information of the product (Note: for each lot of the test kit, this information only needs to be recorded once).
- 3) Check whether the package is intact or not, do not use if it is damaged or leaking.
- 4) Take out the test card from the package and use it within 30 minutes after opening.
- 5) For serum/plasma or whole blood samples, use a pipette to take a certain amount of the sample to be tested and the sample diluent for 1:1 dilution and mix well. Take 75 μL of sample and mix well with the sample diluent 75 μL. Then take 75 μL diluted sample and add to the sample injecting well of the test card directly.
- 6) If the "standard test" mode is selected, the test card should be inserted into the Immunofluorescence Analyzer immediately after the sample is added to the sample injecting well. The analyzer will automatically scan the QR code on the test card and get information. Please select "dilution 1:1" for dilution factor and "serum/plasma" or "whole blood" on the Immunofluorescence Analyzer and press "OK" when first performing the test. For further tests, the analyzer will memorize previous options and the users only need to press "OK" to confirm. Then the device will enter the countdown state and report the result automatically.
- 7) If the "quick test" mode is selected, the timing process should be performed by users. After the sample is added to the sample injecting well, the reaction will take 20 min before the test card can be inserted into the Immunofluorescence Analyzer. At 20 minutes, press the "quick test" mode and insert the test card into the Immunofluorescence Analyzer. Please select "dilution 1:1" for dilution factor and "serum/plasma" or "whole blood" on the Immunofluorescence Analyzer and press "OK" when first performing the test. For further tests, the analyzer will memorize previous options and the users only need to press "OK" to confirm. Under the "quick test" mode, the Immunofluorescence Analyzer will read the product information according to the QR code on the test card and report the test results immediately.
- 8) After the detection is done, the transfer pipette tip, used test card and remaining samples shall be disposed as infectious objects.

【Confidence Interval】

The serum, plasma and whole blood samples of healthy donors (n=226) were tested using the Creatine kinase-MB Test Kit (Quantitative Fluorescent Immuno-chromatography). The results showed that the 95 % percentile of Creatine kinase-MB was 6.4ng/ml.

The above reference values are for reference only. Due to differences in geography, race, sex and age, it is recommended that each laboratory establish its own reference values (reference ranges) according to the actual situation.

【Interpretation of test results】

The results of this reagent are for clinical reference only. The clinical diagnosis and treatment of patients should be considered in combination with their symptoms / signs, medical history, other laboratory examinations and treatment reactions.

【Limitations of test methods】

1. The kit can be used for the determination of serum, plasma or whole blood samples, and the reliability for the concentration determination of Creatine kinase-MB needs further confirmation.
2. Plasma and whole blood samples in this kit are recommended to use EDTA, sodium heparin or sodium citrate as anticoagulants. If other anticoagulants are selected, the effect needs to be further confirmed.
3. The concentrations of Creatine kinase-MB obtained by other methods are not directly comparable with the results of this kit.
4. If the test results of this kit are inconsistent with the clinical evaluations, further examination is required.
5. Operation must be carried out strictly in accordance with operating manual, please operate carefully for getting more accurate results, any modification of the operating manual may affect the results.
6. If false positive results occur, it may be that some non-specific components of the sample that have similar antigenic determinants participate in the recognition reaction; If false negative results occur, it may be that some unknown components block the antigenic determinant and hinder the immune recognition reaction. The detection of Creatine kinase-MB might be unsuccessful if the sample has been taken for too long time or the environmental temperature is causing the degradation of analytes.
7. Other factors may also cause inaccuracy in Creatine kinase-MB test results, including technical problems, inappropriate operation or sample factors, such as wrong sample collection process, preservation and handling.
8. The test results are for reference only and should not be taken as the only indicator for clinical diagnosis.


【Product Performance Requirements】

1. Linear Range within the range of 1 ng/mL to 200 ng/mL, the correlation coefficient R value should be greater than or equal to 0.990.
2. Precision:
 - 2.1 Intra-batch precision: the coefficient of variation (CV) should not be more than 15%.
 - 2.2 Inter-batch precision: the coefficient of variation (CV) shall not be greater than 20%.
3. Accuracy: the relative deviation should not be more than 20%.
4. The limit of detection: The limit of detection of this kit is no higher than 1ng/mL.

【Precautions】

1. This kit is only used for in vitro diagnosis and should not be used for other purposes.
2. This kit is only for use by laboratory medicine professionals.
3. This product is disposable and cannot be reused.
4. Do not use the contaminated or spoiled samples or expired products for testing.
5. Fresh samples are recommended. Hyperlipidemia, jaundice and hemolysis are not recommended.
6. Do not use for detection if the package is not intact.
7. When using the test kit for detection, the environment temperature should not be too high or too low. The test card or test sample stored under low temperature shall recover to room temperature (20 ~ 30 °C) before testing. The test kit should be used within 30 minutes after opening. Please keep dry, and do not exposure the test kit to strong volatile organic reagent (chloroform, toluene and etc.).
8. During the test, the test card and Immunofluorescence Analyzer should avoid strong electromagnetic environment.
9. During the test, do not touch the white film area in the middle of the test card and avoid of touching other contaminants.
10. Do not insert the test card with extra liquids on surface into the analyzer to avoid contaminating or damage of the device.
11. All wastes shall be disposed in accordance with infectious materials.
12. The test result of this kit is not the only indicator for disease. Judgement should be made by doctors while considering all relevant indications and clinical symptoms.
13. It is recommended to wear gloves without talcum powder, because the powder may lead to wrong results.
14. If you have any questions or suggestions about this reagent, please contact the manufacturer.
15. The longest transportation period for the test kit must not exceed 12 days under 2-30 °C.

Index or Symbols

| | | | | | |
|---|-------------------------------------|---|---------------|---|------------------|
|  | Attention, see instructions for use |  | Tests per Kit |  | Manufacturer |
|  | For in vitro |  | Use by |  | Do not reuse |
|  | Store between 2 - 30 °C |  | Lot Number |  | Catalogue Number |

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