## β-human Chorionic Gonadotropin Test Kit (Quantitative Fluorescent Immuno-chromatography)

### 【Product Name】

Generic Name: β-human Chorionic Gonadotropin Test Kit (Quantitative Fluorescent Immuno-chromatography)

### 【Package Spec】

25 PCS/box

### 【Intended Use】

β-Human Chorionic Gonadotropin (β-hCG) Test Kit is used for quantitative detection of β-hCG in female plasma and whole blood samples in vitro.

β-hCG is a glycoprotein secreted by trophoblastic cells of the placenta, which is composed of α and β subunits. The α-subunit is shared by anterior pituitary hormone. The β-subunit is specific for β-hCG. The main function of β-hCG is to stimulate the corpus luteum for continuous secretion of estrogen and progesterone, which in turn promoting formation of decidua in the uterus and maturity of placenta. It is believed β-hCG is produced by trophoblast transitional cells and syncytial cells. β-hCG proliferates rapidly in the first 8 weeks of pregnancy to maintain pregnancy. After around 8 weeks of pregnancy, β-hCG gradually decreases until it reaches relative stability at about 20 weeks.

### 【Principle of Testing】

The Test Kit adopts fluorescence labeling technology. Mouse anti-human β-hCG monoclonal antibody was coated on nitrocellulose membrane. Fluorescence labeled mouse anti-human β-hCG monoclonal antibody was distributed in glass cellulose membrane. If the blood samples contain β-hCG, the analyte can be captured by fluorescence labeled mouse anti-human β-hCG monoclonal antibody. The conjugate can move to the test line area on the nitrocellulose membrane by capillary action and be captured by the corresponding antibody contained therein to form a fluorescent test line. Within a certain range, the concentration of β-hCG in the sample has a positive correlation with the fluorescence intensity of the test line. The concentrations of β-hCG in the sample can be calculated and determined by the immunofluorescence analyzer.

### 【Main components】

The main components of the kit include Test Card, Calibration QR Code (on the box), Diluent and Instruction Manual.

Details are as follows:

1. The test card is composed of a test strip and a plastic clip. The main components of the strip are:

- 1) Mouse anti-human β-hCG antibody (immobilized on the "Test Line" area of nitrocellulose membrane);
- 2) Fluorescent labeled mouse anti-human β-hCG monoclonal antibody (immobilized on the glass cellulose membrane);
- 3) Test strip back supporting material (Polyvinyl Chloride, PVC).

2. Calibration QR code: QR Code located on the box is used to carry the product parameters. Please scan the QR code on the box with Immunofluorescence Analyzer for the input of test information.

3. Diluent

Sample diluent solution: Contains Phosphate buffer and Proclin-300.

Note: Please do not exchange the components in different batches of kits.

### 【Storage conditions and Period of Validity】

The period of validity of the test kit is 24 months while the package is sealed and stored at 2-30 °C. Once the package is opened, please place the test kit at 20-30°C. The opened test kit should be used within 30 min if the humidity is no less than 20% and below 80%RH.

### 【Applicable instruments】

The Immunofluorescence Analyzer manufactured by Deniz Biyokimya Urunleri A.S. R<sub>2</sub>D<sub>2</sub>-V<sub>2</sub>

### 【The requirements of sample】

This kit is suitable for serum/plasma/whole blood samples.

Serum/plasma collection: the serum or plasma should be separated for analysis immediately after blood collection to avoid hemolysis. Use within 24 hours if the sample is stored at room temperature; use within 72 hours if the sample is stored at 2-8 °C; use within 3 months if the sample is stored at -20 °C with one cycle of freeze and thaw.

The collection of whole blood: collect the blood using commercially available blood collection tubes containing EDTA, sodium heparin or sodium citrate anticoagulants. Once the blood is collected, gently invert tube several times to avoid any hemolysis. If the sample cannot be tested within 24 hours, it shall be stored at 2-8 °C. Samples being stored more than 3 days should not be used for examination.

Whole blood/plasma is recommended to be collected using the blood collection tubes containing anticoagulant of EDTA (concentration 3.7-5.4 μmol/mL), sodium heparin (concentration of 12 IU-30 IU/mL of blood) or sodium citrate (concentration 105-109 mmol/L). To separate plasma, it is recommended to centrifuge the blood sample with 500g (3000 rpm) for 5 minutes.

Fresh samples are recommended. Blood samples of hyperlipidemia, jaundice and hemolysis are not recommended.

The sample must reach the equilibrium status at room temperature (20 ~ 30 °C) before use.

### 【Test Method】

Please read the user manual for test kit and Immunofluorescence Analyzer carefully before use, following the requirements and steps strictly. All the test shall be operated under room temperature (20 ~ 30°C) and the relative humidity should be no less than 20%RH and below 80%RH.

The signal ratio (T/C) of Test Line (T) and Quality Control Line (C) reading from the immunofluorescent analyzer is converted to analyte concentration using the standard curve information carried by QR code on the box.

The Quality Control Line is the quality control standard of the test kit. If the Quality Control Line has no signal or abnormal signal (exceed the quality control scope of the analyzer), the test result will be invalid, and it is recommended to repeat the test.

Immunofluorescence Analyzer has two detection modes: "Quick Test" and "Standard Test". Please select the right detection mode according to user manual of Immunofluorescence Analyzer before detection.

The detail operation step is as below:

1) Allow both the test kits and samples balance to room temperature (20 ~ 30 °C) before performing the test. Confirm that the Immunofluorescence Analyzer is connected to powered supply and turned on before use.

2) Before using each lot of the test kit, please select "quick test" mode of the Immunofluorescence Analyzer and read QR code on the box to record the lot information of the product (Note: for each lot of the test kit, this information only needs to be recorded once).

3) Check whether the package is intact or not, do not use if it is damaged or leaking.

4) Take out the test card from the package and use it within 30 minutes after opening.

5) For serum/plasma or whole blood samples, use a pipette to take a certain amount of the sample to be tested and the sample diluent for 1:1 dilution and mix well (E.g. take 75 μL of sample and mix well with the sample diluent 75 μL). Then take 75 μL diluted sample and add to the sample injecting well of the test card directly.

6) When β-hCG exceeds 1000 mIU/mL, should use a pipette to take a certain amount of the sample to be tested and the sample diluent for 1:20 dilution and mix well (E.g. take 20 μL of sample and mix well with the sample diluent 400 μL). Then take 60 μL diluted sample and add to the sample injecting well of the test card directly.

7) If the "Standard Test" mode is selected, the test card should be inserted into the Immunofluorescence Analyzer immediately after the sample is added to the sample injecting well. The analyzer will automatically scan the QR code on the test card and get information. Please select "dilution 1:2" or "dilution 1:20" for dilution factor and "serum/plasma" or "whole blood" on the Immunofluorescence Analyzer and press "OK" when first performing the test. For further tests, the analyzer will memorize previous options and the users only need to press "OK" to confirm. Then the device will enter the countdown state and report the result automatically.

8) If the "Quick Test" mode is selected, the timing process should be performed by users. After the sample is added to the sample injecting well, the reaction will take 15 min before the test card can be inserted into the Immunofluorescence Analyzer and press "OK" when first performing the test. For further tests, the analyzer will memorize previous options and the users only need to press "OK" to confirm. Under the "quick test" mode, the Immunofluorescence Analyzer will read the product information according to the QR code on the test card and report the test results immediately.

9) After the detection is done, the transfer pipette tip, used test card and remaining samples shall be disposed as infectious objects.

### 【Confidence Interval】

The serum, plasma and whole blood samples of non-pregnant women (n=180) were tested using β-human Chorionic Gonadotropin Test Kit (Quantitative Fluorescent Immuno-chromatography). The results showed that the 95 % percentile of β-hCG was <5mIU/mL, the 99 % percentile of β-hCG was <25mIU/mL.

The above reference values are for reference only. Due to differences in geography, race, gender and age, it is recommended that each laboratory establish its own reference interval based on actual conditions.

### 【Interpretation of test results】

The results of this reagent are for clinical reference only. The clinical diagnosis and treatment of patients should be considered in combination with their symptoms / signs, medical history, other laboratory examinations and treatment reactions.

### 【Limitations of test methods】

1. This kit can be used for the determination of serum, plasma or whole blood samples, and the reliability of other humoral samples needs to be determined.

2. Plasma and whole blood samples in this kit are recommended to use EDTA and sodium heparin as anticoagulants. If other anticoagulants are selected, the effect needs to be further confirmed.

3. The concentrations of β-hCG obtained by other methods are not directly comparable with the results of this kit.

4. If the test results of this kit do not match the clinical evaluation, it is recommended to carry out further examination.

5. Operation must be carried out strictly in accordance with operating manual, please operate carefully for getting more accurate results, any modification of the operating manual may affect the results.

6. If false positive results occur, it may be caused by some non-specific components of the sample that have similar antigenic determinants participate in the recognition reaction. In case of false negative results, some unknown components might the antigen-determining cluster and disturb the immune recognition reaction. False negative results may be also caused by deterioration of β-hCG with long waiting time and temperature increase, thus cannot be recognized by the antibody.

7. Other factors may also cause errors in the test results of β-hCG, including technical problems, inappropriate operation or sample factors, such as wrong sample collection process, preservation and handling.

8. The test results are for reference only and are not the only basis for clinical diagnosis.

### 【Product Performance Requirements】

1. Linear range:

1.1 Dilution 1:2: 5 to 1000 mIU/mL, the correlation coefficient r value should be ≥0.990.

1.2 Dilution 1:20: 100 to 20000 mIU/mL, the correlation coefficient r value should be ≥0.990.

2. Precision:

2.1 Repeatability: the coefficient of variation (CV) should not be more than 15%.

2.2 Lot-to-Lot variation: the coefficient of variation (CV) shall not be greater than 20%.

3. Accuracy: the relative deviation should not be more than 20%.

4. The limit of detection: The limit of detection of this kit is no higher than 2.5mIU/mL.

### 【Precautions】

1. This kit is only used for in vitro diagnosis and should not be used for other purposes.

2. The kit is only for use by laboratory medicine professionals.

3. This product is disposable and cannot be reused.

4. Do not use the contaminated or spoiled samples or expired products for testing.

5. Fresh samples are recommended. Hyperlipidemia, jaundice and hemolysis are not recommended.

6. Do not use for detection if the package is not intact.

7. When using the test kit for detection, the environment temperature should not be too high or too low. The test card and test sample stored under low temperature shall recover to room temperature (20 ~ 30 °C) before testing. The test kit should be used within 30 minutes after opening. Please keep dry, and do not exposure the test kit to strong volatile organic reagent (chloroform, toluene and etc.).

8. During the test, the test card and Immunofluorescence Analyzer should avoid strong electromagnetic environment.

9. During the test, do not touch the white film area in the middle of the test card and avoid of touching other contaminants.

10. Do not insert the test card with extra liquids on surface into the analyzer to avoid contaminating or damage of the device.

11. All wastes shall be disposed in accordance with infectious materials.

12. The test result of this kit is not the only indicator for disease. Judgement should be made by doctors while considering all relevant indications and clinical symptoms.

13. It is recommended to wear gloves without talcum powder, because the powder may lead to wrong results.

14. If you have any questions or suggestions about this reagent, please contact the manufacturer.

15. The longest transportation period for the test kit must not exceed 12 days under 2-30 °C.

### Index or Symbols

	Attention, see instructions for use		Tests per Kit		Manufacturer
	For in vitro		Use by		Do not reuse
	Store between 2 - 30°C		Lot Number		Catalogue Number

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